

Air Regulator / Filter

Use and mounting directions



Use this regulator / filter to connect the Scale–Aire, Scale–Aire Mini, and the Drill–Aire Series, and other dental air units, to a compressor.

Use the supplied mounting hardware to attach the regulator to drywall.

The unit must be mounted on a wall, with the clear reservoir down to function properly.

Caution: Mounting this regulator incorrectly can damage your equipment and void your warranty.

The reservoir is clear to allow for visual inspection of accumulated water. A spigot is located at the bottom of the reservoir for convenient draining. Drain daily.

Use the supplied male quick connect / hose to connect to your dental air unit, connect the other hose to your compressor.

This unit is pre-calibrated to 75 PSI prior to shipping. Exceeding 80 PSI can cause serious injury, damage the equipment and void the warranty

CAUTION: This regulator is not suitable for connecting to a nitrogen cylinder.

Call Engler Engineering Corp. at 800-445-8581 for information on connecting a Scale – Aire Mini to a nitrogen tank (cylinder)



Web: www.englerusa.com / Email: info@englerusa.com
www.dynaxusa.com / www.ez-stretcher.com / www.ez-dentaltable.com /
www.25KSeries.com for brochures, manuals and "how to"
pages please visit www.engler411.com

Engler Engineering Corporation
1099 East 47th Street - Hialeah, Florida 33013 USA
Tel: (800) 445-8581 / (305) 688-8581 / Fax: (305) 685-7671



Proudly made
in the U.S.A.