


# MAXITHERM<sup>®</sup>

## Heavy Duty Reusable Vinyl Circulator Pads


Instructions for connecting the Maxitherm<sup>®</sup> to your Gaymar<sup>®</sup> and most other circulators.

Step 1: Remove the old pad from your circulator hose. Remove the Gaymar<sup>®</sup> connectors (both male & Female) from the hose by immersing the connectors into hot water. It may be necessary to cut the connectors off of the hose (save the connectors for future use).



Step 2: The Maxitherm<sup>®</sup> Pad has 2 male connectors intended for the Microtemp<sup>®</sup> circulator. It will be necessary to remove these connectors. Each Maxitherm<sup>®</sup> contains a small package with 4 extra connectors. The two female parts can be inserted into the Gaymar<sup>®</sup> hose while the two male parts will replace the connectors on the Maxitherm<sup>®</sup> pad.

Sm	12" X 18"
Med	22" X 30"
Lg	24" X 60"



Step 3: At this time, your new Maxitherm<sup>®</sup> circulator pad will connect directly to your Gaymar<sup>®</sup> hose and circulator.

\* Gaymar<sup>®</sup> is the registered trademark of GAYMAR INDUSTRIES INC.