


Air Regulator / Filter

Mounting and Operating Instructions



Use this regulator / filter to connect a compressor to the Engler Scale-Aire Mini / Drill-Aire Series.

The unit must be mounted vertically, with the clear reservoir pointed down to function properly.

Caution: Mounting this regulator incorrectly can damage your equipment and void your warranty.

The reservoir is clear to allow for visual inspection of accumulated water. A spigot (drain valve) is located at the bottom of the clear reservoir for convenient draining. Drain daily.

Use the supplied male quick connect / hose to connect to your dental air unit, connect the other hose to your compressor. Additional optional hardware may be required to connect the compressor. See video link below.

Maximum **INPUT** pressure: 175 PSI (12 BAR). Maximum temperature: 125F (80C)

This unit is pre-calibrated to 75 PSI **OUTPUT** prior to shipping. Exceeding 80 PSI can cause serious injury, damage the equipment (handpieces) and void the warranty.

CAUTION: This regulator is not suitable for connecting to a nitrogen cylinder.

For video instructions please visit; <https://engler411.com/youtube/videos08.html>, then click on video #4



Web: www.englerusa.com / Email: info@englerusa.com
www.dynaxusa.com / www.englerO2.com
www.ez-breathe.com / www.25KSeries.com
for brochures, manuals and "how to"
pages please visit www.engler411.com



Proudly made
in the U.S.A.

Engler Engineering Corporation
1099 East 47th Street - Hialeah, Florida 33013 USA
Tel: (800) 445-8581 / (305) 688-8581 / Fax: (305) 685-7671